Wise Wave IT Solutions

LIMS

Wise Wave IT Solutions offers a range of solutions designed to streamline your business operations and enhance efficiency.



Wise Wave IT Solutions' Specialization



Custom Software Solutions

At Wise Wave IT Solutions we specialize in crafting Customized solutions that align perfectly with your business objective and unique challenges.

By recognizing that every organization has its own identity, we embrace a collaborative approach to design and develop innovative software that derives your success and transforms your vision into reality.

Key Features

Agile Development

Flexible and iterative approach ensures raped development and quick to changes

User Centric Design

We optimize user experience, crafting innovative interface that are easy to navigate.

Scalable Solutions

At Our software is designed to grow with your business, adapting to changing needs and seamlessly integrating new features.

Laboratory Information Management System(LIMS)

* A Standard & Advanced LIMS Solution for Analytical Laboratories

• ISO 17025 Compliance through Our LIMS

Our Laboratory Information Management System (LIMS) seamlessly integrates all the essential requirements of ISO 17025, guaranteeing complete adherence to globally recognized standards for testing and calibration laboratories.

NABL Accreditation Support through Our LIMS

With robust features tailored to meet NABL accreditation standards, our LIMS ensures your laboratory consistently adheres to rigorous guidelines, enhancing accuracy, traceability, and reliability across all testing and calibration processes..

Our Laboratory Information Management System (LIMS) is designed to enhance efficiency, ensure compliance, and optimize workflows across all laboratory functions.

Advantages of our LIMS

Improved Efficiency

LIMS automates and streamlines laboratory processes, reducing manual tasks and human error, leading to faster turnaround times for tests and reports.

Increased Productivity

By automating routine tasks such as sample tracking and result reporting, LIMS frees up time for laboratory staff to focus on high-priority tasks, increasing overall productivity.

Cost Savings

LIMS helps reduce operational costs by minimizing errors, reducing paper-based processes, and improving resource utilization, leading to better overall cost management.

Scalability

As your laboratory grows, LIMS can easily scale to accommodate additional samples, tests, and users, ensuring continued efficiency without significant extra overhead.

• Enhanced Data Management

It centralizes data storage, making it easy to manage, access, and retrieve test results, improving accuracy and consistency across the laboratory.

Compliance and Traceability

LIMS ensures compliance with industry standards like ISO 17025 and NABL by maintaining detailed records, audit trails, and traceability for every sample tested.

Accurate Reporting


LIMS provides customizable reporting tools, allowing laboratories to generate accurate, ondemand reports that meet regulatory requirements and client expectations.


Enhanced Security

With built-in security features, LIMS protects sensitive data through encryption, access control, and user authentication, ensuring data integrity and confidentiality.

Features of LIMS







Manage Masters Masters like
Department,
Employee, customer,
Test are needs to
feed one time in
system, value of
these masters are
used in another
modules

Sample Registration SRF No./ work Order No., Sample No. Is generated on filling up the Sample Request / work Order form.

Sample Review Sample Review by QM/TM. QM/TM can Accept or Reject sample based on sample review criteria.

Job Allocation QM/ TM can allocate sample wise test parameter to department or analyst.

Test Perform Analyst performs test and feeds observation, system will do all the calculations and validate result as per the master criteria.



QM will review the test result feed by analyst and approve/return based on observation, On approval ULR No will be generated for sample having NABL test parameters.



Authorized user can take print and send test reports to client.



Account department user can generate quotation/invoice/generate payment receipt

Inventory Manageme nt

Inventory department can manage inventory by inward, outward and stock reports

Equipment Calibration

System manages equipment wise calibration history and notify respective users for next calibration



System manages all the NABL audit required documents




System generates various reports including customized reports



System manages customer feedback

LIMS Training and Support



Technical details our LIMS

- Front End
 - Microsoft Blazor
 - Net Core 8.0
 - .Net Core Web API
- Hosting Server Configuration
 - Windows Server 2019 or Above
 - ❖ IIS Server Version = 10.0
 - Net Core 8.0
 - VM 16 Core, 32 GB RAM
 - ❖ Minimum 150 GB HDD
 - ❖ MySQL version 8.4/ MSSQL server version 2022

- Back End
 - ❖ MySQL/MS SQL Server

- Client System Configuration
 - Windows 10 or Above
 - ❖ *VM 4 Core, 8 GB RAM*
 - Browser Support: Chrome, Edge, Firefox

Online /Offline Hosting of our LIMS

Online Cloud-base Hosting

- Remote Accessibility
- Scalability
- Customization
- Automatic Updates
- Data Backup and Disaster Recovery
- Collaboration and Integration of multiple labs
- Security

Online SAAS

- Remote Accessibility
- Automatic Updates & Maintenance
- Cost Efficiency
- Flexible Payment Options
- No Upfront Capital Expenses
- Security and Data Privacy
- Lab and user wise subscription

Offline Hosting(On-premise Hosting)

- Data Control and Security
- No Dependency on Internet Connectivity
- Customization
- One-time Cost
- Faster Performance
- No Vendor Lock-In



Connect with Us

wisewaveitsolutions.com

Email: info@wisewaveitsolutions.com

Call on: +91-7073203127

Visit at : G Floor, 40/C, Parivahan Nagar, Khatipura, Jaipur, Rajasthan-302012, India